

Anti-PAFAH1B3 Antibody

Catalog Number: A07650

About PAFAH1B3

PAFAH1B3, also named as PAFAHG, belongs to the 'GDSSL' lipolytic enzyme family and Platelet-activating factor acetylhydrolase IB beta/gamma subunits subfamily. It inactivates paf by removing the acetyl group at the sn-2 position. This is a catalytic subunit. PAFAH1B3 plays an important role during the development of brain. The antibody is specific to PAFAH1B3.

Overview

Product Name	Anti-PAFAH1B3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PAFAH1B3 Antibody catalog # A07650. Tested in WB, IHC, ICC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q15102

Technical Details

Immunogen	Peptide
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PAFAH1B3 Antibody

For Research Use Only. Not for use in diagnostic procedures.