

## Anti-PGM2 Antibody Picoband® Cy3 Conjugated

Catalog Number: A07690-1-Cy3

### About PGM2

Phosphoglucomutase-2 is an enzyme that in humans is encoded by the PGM2 gene.

### Overview

Product Name	Anti-PGM2 Antibody Picoband® Cy3 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q96G03

### Technical Details

Immunogen	E.coli-derived human PGM2 recombinant protein (Position: R72-D520). Human PGM2 shares 92.7% amino acid (aa) sequence identity with mouse PGM2.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-PGM2 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.