

Anti-RHBDF2 Antibody Picoband® Biotin Conjugated

Catalog Number: A07693-1-Biotin

About RHBDF2

Rhomboid family member 2 is a protein that in humans is encoded by the RHBDF2 gene. Predicted to enable protein transporter activity. Predicted to be involved in negative regulation of protein secretion and regulation of epidermal growth factor receptor signaling pathway. Located in plasma membrane. Implicated in palmoplantar keratoderma-esophageal carcinoma syndrome.

Overview

Product Name	Anti-RHBDF2 Antibody Picoband® Biotin Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q6PJF5

Technical Details

Immunogen	E.coli-derived human RHBDF2 recombinant protein (Position: H145-D529).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RHBDF2 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.