

Anti-RHBDF2 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A07693-1-Fluoro647

About RHBDF2

Rhomboid family member 2 is a protein that in humans is encoded by the RHBDF2 gene. Predicted to enable protein transporter activity. Predicted to be involved in negative regulation of protein secretion and regulation of epidermal growth factor receptor signaling pathway. Located in plasma membrane. Implicated in palmoplantar keratoderma-esophageal carcinoma syndrome.

Overview

Product Name	Anti-RHBDF2 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6PJF5

Technical Details

Immunogen	E.coli-derived human RHBDF2 recombinant protein (Position: H145-D529).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RHBDF2 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.