

## Anti-Sarcalumenin/Srl Antibody Picoband® HRP Conjugated

Catalog Number: A07699-4-HRP

### About Srl

Sarcalumenin is a protein that in humans is encoded by the SRL gene. Predicted to enable GTP binding activity. Acts upstream of or within response to muscle activity involved in regulation of muscle adaptation and store-operated calcium entry. Predicted to be located in membrane and sarcoplasmic reticulum. Predicted to be active in cytoplasm; intracellular membrane-bounded organelle; and plasma membrane. Is expressed in brain ventricular layer; diaphragm; limb; skeletal musculature; and tail mesenchyme. Orthologous to human SRL (sarcalumenin).

### Overview

Product Name	Anti-Sarcalumenin/Srl Antibody Picoband® HRP Conjugated
Reactive Species	Mouse
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q7TQ48

### Technical Details

Immunogen	E.coli-derived mouse Sarcalumenin/Srl recombinant protein (Position: R119-M422).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-Sarcalumenin/Srl Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.