

Anti-RFC4 Antibody Picoband® HRP Conjugated

Catalog Number: A07702-1-HRP

About RFC4

Replication factor C subunit 4 is a protein that in humans is encoded by the RFC4 gene. The elongation of primed DNA templates by DNA polymerase delta and DNA polymerase epsilon requires the accessory proteins proliferating cell nuclear antigen (PCNA) and replication factor C (RFC). RFC, also named activator 1, is a protein complex consisting of five distinct subunits of 140, 40, 38, 37, and 36 kD. This gene encodes the 37 kD subunit. This subunit forms a core complex with the 36 and 40 kDa subunits. The core complex possesses DNA-dependent ATPase activity, which was found to be stimulated by PCNA in an in vitro system. Alternatively spliced transcript variants encoding the same protein have been reported.

Overview

| | |
|----------------------|---|
| Product Name | Anti-RFC4 Antibody Picoband® HRP Conjugated |
| Reactive Species | Human |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. |
| Host | Rabbit |
| Uniprot ID | P35249 |

Technical Details

| | |
|---------------------|--|
| Immunogen | E.coli-derived human RFC4 recombinant protein (Position: E42-Q323). |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | HRP |
| Suggested Dilutions | The intended application should be selected according to the customer's experimental requirements. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RFC4 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.