

Anti-ATG101 Antibody Picoband® Biotin Conjugated

Catalog Number: A07713-2-Biotin

About ATG101

Enables identical protein binding activity. Involved in autophagosome assembly and negative regulation of cell population proliferation. Located in phagophore assembly site. Part of Atg1/ULK1 kinase complex. Biomarker of cholangiocarcinoma; colon cancer; and hepatocellular carcinoma.

Overview

Product Name	Anti-ATG101 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9BSB4

Technical Details

Immunogen	E.coli-derived human ATG101 recombinant protein (Position: Q53-A199).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ATG101 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.