

## Anti-Zinc transporter SLC39A7 SLC39A7 Antibody

Catalog Number: A07719-3


### Overview

Product Name	Anti-Zinc transporter SLC39A7 SLC39A7 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Zinc transporter SLC39A7 SLC39A7 Antibody catalog # A07719-3. Tested in WB,ICC applications. This antibody reacts with Human,Mouse,Rat.
Application	ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q92504

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 179-224 of Human SLC39A7.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 ICC: 1:50-1:200

## Anti-Zinc transporter SLC39A7 SLC39A7 Antibody (A07719-3) Images



Western blot (WB) analysis of SLC39A7 polyclonal antibody at 1:500 dilution Lane1:Hela cell lysate Lane2:sp2/0 cell lysate Lane3:H9C2 cell lysate

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zinc transporter SLC39A7 SLC39A7 Antibody

For Research Use Only. Not for use in diagnostic procedures.