

## Anti-TEAD3 Antibody Picoband® Biotin Conjugated

Catalog Number: A07748-2-Biotin

### About TEAD3

Transcriptional enhancer factor TEF-5 is a protein that in humans is encoded by the TEAD3 gene. This gene product is a member of the transcriptional enhancer factor (TEF) family of transcription factors, which contain the TEA/ATTS DNA-binding domain. It is predominantly expressed in the placenta and is involved in the transactivation of the chorionic somatomammotropin-B gene enhancer. Translation of this protein is initiated at a non-AUG (AUA) start codon.

### Overview

Product Name	Anti-TEAD3 Antibody Picoband® Biotin Conjugated
Reactive Species	Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q99594

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human TEAD3, identical to the related mouse sequences.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TEAD3 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.