

Anti-RNMT Antibody Picoband® Biotin Conjugated

Catalog Number: A07752-1-Biotin

About RNMT

Enables RNA binding activity and mRNA (guanine-N7-)-methyltransferase activity. Involved in 7-methylguanosine mRNA capping. Located in fibrillar center and nucleoplasm. Part of mRNA cap binding activity complex; mRNA cap methyltransferase complex; and receptor complex.

Overview

Product Name	Anti-RNMT Antibody Picoband® Biotin Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O43148

Technical Details

Immunogen	E.coli-derived human RNMT recombinant protein (Position: R197-E380).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-RNMT Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.