

Anti-PPP1R3A Antibody

Catalog Number: A07755


Overview

| | |
|----------------------|--|
| Product Name | Anti-PPP1R3A Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-PPP1R3A Antibody (Catalog# A07755). Tested in WB, ELISA applications. This antibody reacts with Human, Rat, Mouse. |
| Application | ELISA, WB |
| Clonality | Polyclonal |
| Formulation | PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage Instructions | Store at -20°C for one year. For short-term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | Q16821 |

Technical Details

| | |
|----------------------------|--|
| Immunogen | A peptide derived from human PPP1R3A. Immunogen sequence location: 647 - 696 |
| Predicted Reactive Species | Bovine, Canine, Chicken, Primate, Sheep, Xenopus, Zebrafish |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure. |
| Purification | The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen. |
| Suggested Dilutions | <p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows:</p> <p>Boster Bio's internal QC testing used:</p> <p>Western Blot, 1:500 - 1:2000</p> <p>ELISA, 1:20000</p> <p>Not yet tested in other applications.</p> |

Anti-PPP1R3A Antibody (A07755) Images



Western blot analysis of lysates from MCF-7 and COLO cells, using PPP1R3A Antibody. The lane on the right is blocked with the synthesized peptide.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PPP1R3A Antibody