

Anti-SLC32A1 Antibody

Catalog Number: A07775


Overview

Product Name	Anti-SLC32A1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-SLC32A1 Antibody catalog # A07775. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9H598

Technical Details

Immunogen	Recombinant fusion protein of human SLC32A1(Q9H598).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000-1:2000

Anti-SLC32A1 Antibody (A07775) Images



Western blot analysis of extracts of various cell lines, using SLC32A1 antibody at 1:400 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.

1 Publications Citing This Product

1. PubMed ID: 10.1038/s41598-019-51880-9, A unique olfactory bulb microcircuit driven by neurons expressing the precursor to glucagon-like peptide 1

Visit bosterbio.com/anti-slc32a1-antibody-a07775-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SLC32A1 Antibody

For Research Use Only. Not for use in diagnostic procedures.