

Anti-Oip5 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A07819-1-Fluoro488

About Oip5

Predicted to enable identical protein binding activity. Predicted to be involved in CENP-A containing chromatin assembly and chromosome segregation. Predicted to act upstream of or within cell division. Located in chromocenter and nucleus. Orthologous to human OIP5 (Opa interacting protein 5).

Overview

Product Name	Anti-Oip5 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	A2AQ14

Technical Details

Immunogen	E.coli-derived mouse Oip5 recombinant protein (Position: M1-K220).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Oip5 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.