

Anti-Oip5 Antibody Picoband® HRP Conjugated

Catalog Number: A07819-1-HRP

About Oip5

Predicted to enable identical protein binding activity. Predicted to be involved in CENP-A containing chromatin assembly and chromosome segregation. Predicted to act upstream of or within cell division. Located in chromocenter and nucleus. Orthologous to human OIP5 (Opa interacting protein 5).

Overview

Product Name	Anti-Oip5 Antibody Picoband® HRP Conjugated
Reactive Species	Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	A2AQ14

Technical Details

Immunogen	E.coli-derived mouse Oip5 recombinant protein (Position: M1-K220).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Oip5 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.