

Anti-ELP2 Antibody Picoband® HRP Conjugated

Catalog Number: A07823-2-HRP

About ELP2

The protein encoded by this gene is a core subunit of the elongator complex, a histone acetyltransferase complex that associates with RNA polymerase II. In addition to histone acetylation, the encoded protein effects transcriptional elongation and may help remodel chromatin.

Overview

Product Name	Anti-ELP2 Antibody Picoband® HRP Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q6IA86

Technical Details

Immunogen	E.coli-derived human ELP2 recombinant protein (Position: D84-L826). Human ELP2 shares 77.9% and 79.8% amino acid (aa) sequence identity with mouse and rat ELP2, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ELP2 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.