

Anti-DHRS3 Antibody

Catalog Number: A07824

About DHRS3

Predicted to enable NAD-retinol dehydrogenase activity. Predicted to be involved in regulation of retinoic acid receptor signaling pathway and retinoid metabolic process. Predicted to act upstream of or within several processes, including animal organ morphogenesis; negative regulation of retinoic acid receptor signaling pathway; and regulation of ossification. Predicted to be located in endoplasmic reticulum membrane and photoreceptor outer segment membrane. Predicted to be active in lipid droplet.

Overview

Product Name	Anti-DHRS3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-DHRS3 Antibody catalog # A07824. Tested in WB, ICC, IF, IP, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IP, IF, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O75911

Technical Details

Immunogen	E.coli-derived human DHRS3 recombinant protein (Position: R5-T302).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunocytochemistry/Immunofluorescence, 1:50-400 Immunoprecipitation, 1:50 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DHRS3 Antibody

For Research Use Only. Not for use in diagnostic procedures.