

## Anti-OTUD4 Antibody Picoband® Biotin Conjugated

Catalog Number: A07830-2-Biotin

### About OTUD4

Alternatively spliced transcript variants have been found for this gene. The smaller protein isoform encoded by the shorter transcript variant is found only in HIV-1 infected cells.

### Overview

Product Name	Anti-OTUD4 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q01804

### Technical Details

Immunogen	E.coli-derived human OTUD4 recombinant protein (Position: E52-T1114). Human OTUD4 shares 87% amino acid (aa) sequence identity with mouse OTUD4.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-OTUD4 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.