

Anti-NaDC-1/SLC13A2 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A07831-1-Fluoro647

About SLC13A2

Solute carrier family 13 member 2 is a protein that is encoded in humans by the SLC13A2 gene. The protein encoded by this gene is a sodium-coupled citrate transporter that is regulated by the chaperone activity of cyclophilin b. The encoded protein may play a role in the formation of kidney stones.

Overview

Product Name	Anti-NaDC-1/SLC13A2 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q13183

Technical Details

Immunogen	E.coli-derived human NaDC-1/SLC13A2 recombinant protein (Position: A152-M417).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NaDC-1/SLC13A2 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.