

Anti-COQ4 Antibody

Catalog Number: A07857

About COQ4

This gene encodes a component of the coenzyme Q biosynthesis pathway. Coenzyme Q, an essential component of the electron transport chain, shuttles electrons between complexes I or II to complex III of the mitochondrial transport chain. This protein appears to play a structural role in stabilizing a complex that contains most of the coenzyme Q biosynthesis enzymes. Mutations in this gene are associated with mitochondrial disorders linked to coenzyme Q deficiency. Alternative splicing results in multiple transcript variants.

Overview

Product Name	Anti-COQ4 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-COQ4 Antibody catalog # A07857. Tested in WB, IHC, IP applications. This antibody reacts with Human, Mouse, Rat.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9Y3A0

Technical Details

Immunogen	E.coli-derived human COQ4 recombinant protein (Position: T3-Q262).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 ImmunoPrecipitation, 1:250-300

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-COQ4 Antibody

For Research Use Only. Not for use in diagnostic procedures.