

Anti-CNG-2 CNGA2 Antibody

Catalog Number: A07901

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-CNG-2 CNGA2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-CNG-2 CNGA2 Antibody catalog # A07901. Tested in WB, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q16280

Technical Details


Immunogen	The antiserum was produced against synthesized peptide derived from human CNGA2. AA range:391-440
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000

	IF 1:200-1:1000 ELISA 1:40000
--	----------------------------------


Anti-CNG-2 CNGA2 Antibody (A07901) Images



Western Blot analysis of various cells using CNG-2 Polyclonal Antibody



Western Blot analysis of K562 cells using CNG-2 Polyclonal Antibody



Immunofluorescence analysis of A549 cells, using CNGA2 Antibody. The picture on the right is blocked with the synthesized peptide.

Anti-CNG-2 CNGA2 Antibody

For Research Use Only. Not for use in diagnostic procedures.