

Anti-NDUFB8 Antibody

Catalog Number: A07936


Overview

Product Name	Anti-NDUFB8 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-NDUFB8 Antibody catalog # A07936. Tested in WB,IHC,ICC/IF,IP,Flow Cytometry applications. This antibody reacts with Human,Mouse,Rat.
Application	Flow Cytometry, IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O95169

Technical Details

Immunogen	A synthesized peptide derived from human NDUFB8.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Protein A affinity purified.
Suggested Dilutions	WB: 1:500-1:1000 IHC-P: 1:50-1:100 ICC/IF: 1:50-1:100 IP: 1:30 FC: 1:50-1:100

Anti-NDUFB8 Antibody (A07936) Images



Western blot analysis of NDUFB8 in HeLa \square 3T3 \square PC-12 lysates using NDUFB8 antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NDUFB8 Antibody

For Research Use Only. Not for use in diagnostic procedures.