

Anti-Dok-2 (phospho-Y299) Antibody

Catalog Number: A07956Y299


Overview

Product Name	Anti-Dok-2 (phospho-Y299) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Dok-2 (phospho-Y299) Antibody catalog # A07956Y299. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O60496

Technical Details

Immunogen	Synthetic phosphopeptide derived from human Dok-2 around the phosphorylation site of Tyrosine 299.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

Anti-Dok-2 (phospho-Y299) Antibody (A07956Y299) Images



Western blot (WB) analysis of p-Dok-2 (Y299) pAb at 1:500 dilution Lane1:HEK293T whole cell lysate(40ug)
Lane2:HEK293T treated with EGF(100ng/ml,30 minutes) whole cell lysate(40ug) Lane3:H9C2 whole cell lysate(40ug)
Lane4:3T3-L1 whole cell lysate(40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Dok-2 (phospho-Y299) Antibody

For Research Use Only. Not for use in diagnostic procedures.