

Anti-ARHGAP17 (T363) Antibody

Catalog Number: A07968T363


Overview

Product Name	Anti-ARHGAP17 (T363) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ARHGAP17 (T363) Antibody catalog # A07968T363. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q68EM7

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 845-905 of Human RoXaN.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-ARHGAP17 (T363) Antibody (A07968T363) Images



Western blot (WB) analysis of ARHGAP17 (T363) pAb at 1:1000 dilution Lane1:Hela whole cell lysate(30ug) Lane2:HCT116 whole cell lysate(40ug) Lane3:A375 whole cell lysate(40ug) Lane4:The Embryo tissue lysate of Mouse(40ug) Lane5:H9C2 whole cell lysate(40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ARHGAP17 (T363) Antibody

For Research Use Only. Not for use in diagnostic procedures.