

## Anti-ACSS1 Antibody

Catalog Number: A07972


### Overview

Product Name	Anti-ACSS1 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-ACSS1 Antibody catalog # A07972. Tested in WB applications. This antibody reacts with Human, Mouse.
Application	WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9NUB1


### Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human ACSS1. AA range:611-660
Isotype	IgG
Form	Liquid
Concentration	This antibody's concentration is >0.5mg/ml.
Purification	Immunogen affinity purified
Suggested Dilutions	WB, 1:500-1:2000

## Anti-ACSS1 Antibody (A07972) Images



Western Blot (WB) analysis of HT29 cells using ACSS1 Polyclonal antibody.



Western Blot (WB) analysis of specific cells using ACSS1 Polyclonal antibody.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ACSS1 Antibody