

## Anti-Melanoma-associated antigen 6 MAGEA6 Antibody

Catalog Number: A07982


### Overview

|                      |   |
|----------------------|---|
| Product Name         | Anti-Melanoma-associated antigen 6 MAGEA6 Antibody  |
| Reactive Species     | Human, Mouse, Rat   |
| Description          | Boster Bio Anti-Melanoma-associated antigen 6 MAGEA6 Antibody catalog # A07982. Tested in WB applications. This antibody reacts with Human,Mouse,Rat. |
| Application          | WB  |
| Clonality            | Polyclonal  |
| Formulation          | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2  |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.            |
| Host                 | Rabbit  |
| Uniprot ID           | P43360  |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | A synthetic peptide corresponding to residues in Human MAGEA6.  |
| Isotype             | IgG   |
| Form                | Liquid  |
| Concentration       | 1 mg/ml   |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | WB: 1:500-1:1000  |

## Anti-Melanoma-associated antigen 6 MAGEA6 Antibody (A07982) Images



Western blot (WB) analysis of MAGEA6 polyclonal antibody at 1:500 dilution Lane1:MCF-7 cell lysate Lane2:NIH-3T3 cell lysate Lane3:H9C2 cell lysate

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Melanoma-associated antigen 6 MAGEA6 Antibody

For Research Use Only. Not for use in diagnostic procedures.