

Anti-Treslin TICRR Antibody

Catalog Number: A08002

About TICRR

TopBP1-interacting checkpoint and replication regulator, Treslin, chromosome 15 open reading frame 42, C15ORF42


Overview

Product Name	Anti-Treslin TICRR Antibody
Reactive Species	Human
Description	Boster Bio Anti-Treslin TICRR Antibody (Catalog # A08002). Tested in ELISA, WB, IHC-P, IF applications. This antibody reacts with Human.
Application	ELISA, IF, IHC-P, WB
Clonality	Polyclonal
Formulation	TICRR antibody is supplied in PBS containing 0.02% sodium azide.
Storage Instructions	TICRR antibody can be stored at 4°C for three months and -20°C, stable for up to one year.
Host	Rabbit
Uniprot ID	Q7Z2Z1


Technical Details

Immunogen	TICRR antibody was raised against a 19 amino acid peptide near the carboxy terminus of human TICRR. The immunogen is located within amino acids 1840 - 1890 of TICRR.
Predicted Reactive Species	Mouse
Cross Reactivity	TICRR antibody is human specific. At least two isoforms of TICRR are known to exist.
Isotype	IgG
Form	Liquid
Concentration	1 mg/mL
Purification	TICRR antibody is affinity chromatography purified via peptide column.
Suggested Dilutions	TICRR antibody can be used for detection of TICRR by Western blot at 1 - 2 ug/ml. Antibody can also be used for Immunohistochemistry at 5 ug/mL. For Immunofluorescence start at 20 ug/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested. Optimal dilutions for each application should be determined by the researcher.


Anti-Treslin TICRR Antibody (A08002) Images



Western blot analysis of TICRR in human small intestine tissue lysate with TICRR antibody at 1 ug/ml.



Immunohistochemistry of TICRR in human small intestine tissue with TICRR antibody at 5 ug/mL.



Immunofluorescence of TICRR in human small intestine tissue with TICRR antibody at 20 ug/mL.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Treslin TICRR Antibody

For Research Use Only. Not for use in diagnostic procedures.