

Anti-Tryptase Alpha / Beta 1 TPSAB1 Antibody

Catalog Number: A08039-1


Overview

Product Name	Anti-Tryptase Alpha / Beta 1 TPSAB1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Tryptase Alpha / Beta 1 TPSAB1 Antibody catalog # A08039-1. Tested in WB,IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q15661


Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:5,000 IHC: 1:50-1:200

Anti-Tryptase Alpha / Beta 1 TPSAB1 Antibody (A08039-1) Images



Western blot analysis of Mast Cell Tryptase on human skin lysates using anti-Mast Cell Tryptase antibody at 1/1,000 dilution.



Immunohistochemical analysis of paraffin-embedded human tonsil tissue using anti-Mast Cell Tryptase antibody. Counter stained with hematoxylin.

1 Publications Citing This Product

1. PubMed ID: 24735816, Tryptase inhibitor APC 366 prevents hepatic fibrosis by inhibiting collagen synthesis induced by tryptase/protease-activated receptor 2 interactions in hepatic stellate cells

Visit bosterbio.com/anti-tryptase-alpha-beta-1-tpsab1-antibody-a08039-1-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Tryptase Alpha / Beta 1 TPSAB1 Antibody

For Research Use Only. Not for use in diagnostic procedures.