

## Anti-VBP1 Antibody

Catalog Number: A08073-1

### About VBP1

The protein encoded by this gene interacts with the Von Hippel-Lindau protein to form an intracellular complex. The encoded protein functions as a chaperone protein, and may play a role in the transport of the Von Hippel-Lindau protein from the perinuclear granules to the nucleus or cytoplasm. Alternative splicing and the use of alternate transcription start sites results in multiple transcript variants encoding different protein isoforms.

### Overview

Product Name	Anti-VBP1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-VBP1 Antibody catalog # A08073-1. Tested in WB, IHC, ICC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt $-20^{\circ}\text{C}$ as supplied. 6 months $2$ to $8^{\circ}\text{C}$ after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P61758

### Technical Details

Immunogen	E.coli-derived human VBP1 recombinant protein (Position: A3-N195).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400 ELISA, 1:100-1000

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-VBP1 Antibody

For Research Use Only. Not for use in diagnostic procedures.