

Anti-PVRL2 NECTIN2 Antibody

Catalog Number: A08081-1


Overview

Product Name	Anti-PVRL2 NECTIN2 Antibody
Reactive Species	Human
Description	Boster Bio Anti-PVRL2 NECTIN2 Antibody catalog # A08081-1. Tested in WB,IHC applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q92692

Technical Details

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 382-479 of human Nectin 2/CD112 (NP_002847.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:50-1:200

Anti-PVRL2 NECTIN2 Antibody (A08081-1) Images



Western blot analysis of extracts of various cell lines, using Nectin 2/CD112 antibody at 1:1000 dilution.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PVRL2 NECTIN2 Antibody

For Research Use Only. Not for use in diagnostic procedures.