

## Anti-NSFL1 cofactor p47 NSFL1C Antibody

Catalog Number: A08116


### Overview

Product Name	Anti-NSFL1 cofactor p47 NSFL1C Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-NSFL1 cofactor p47 NSFL1C Antibody catalog # A08116. Tested in WB,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UNZ2


### Technical Details

Immunogen	Recombinant full length Human NSFL1C.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IF: 1:50-1:200

## Anti-NSFL1 cofactor p47 NSFL1C Antibody (A08116) Images



WesternBlot (WB) analysis of NSFL1C polyclonal antibody



Immunoprecipitation analysis of 200ug extracts of A-549 cells, using 3 ug NSFL1C antibody.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NSFL1 cofactor p47 NSFL1C Antibody

For Research Use Only. Not for use in diagnostic procedures.