

Anti-Nir2/PITPNM1 Antibody Picoband® Biotin Conjugated

Catalog Number: A08128-2-Biotin

About Nir2/PITPNM1

Membrane-associated phosphatidylinositol transfer protein 1 is a protein that in humans is encoded by the PITPNM1 gene. PITPNM1 belongs to a family of membrane-associated phosphatidylinositol transfer domain-containing proteins that share homology with the Drosophila retinal degeneration B (rdgB) protein.

Overview

Product Name	Anti-Nir2/PITPNM1 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O00562

Technical Details

Immunogen	E.coli-derived human Nir2/PITPNM1 recombinant protein (Position: E418-R831). Human PITPNM1 shares 88.9% and 89.9% amino acid (aa) sequence identity with mouse and rat PITPNM1, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-Nir2/PITPNM1 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.