

## Anti-Nir2/PITPNM1 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A08128-2-Fluoro550

### About Nir2/PITPNM1

Membrane-associated phosphatidylinositol transfer protein 1 is a protein that in humans is encoded by the PITPNM1 gene. PITPNM1 belongs to a family of membrane-associated phosphatidylinositol transfer domain-containing proteins that share homology with the Drosophila retinal degeneration B (rdgB) protein.

### Overview

Product Name	Anti-Nir2/PITPNM1 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	O00562

### Technical Details

Immunogen	E.coli-derived human Nir2/PITPNM1 recombinant protein (Position: E418-R831). Human PITPNM1 shares 88.9% and 89.9% amino acid (aa) sequence identity with mouse and rat PITPNM1, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-Nir2/PITPNM1 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.