

Anti-PGS1 Antibody Picoband® HRP Conjugated

Catalog Number: A08135-HRP

About PGS1

Predicted to enable CDP-diacylglycerol-glycerol-3-phosphate 3-phosphatidyltransferase activity and calcium ion binding activity. Predicted to be involved in cardiolipin biosynthetic process and diacylglycerol metabolic process. Located in endoplasmic reticulum. Is active in mitochondrion.

Overview

Product Name	Anti-PGS1 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q32NB8

Technical Details

Immunogen	E.coli-derived human PGS1 recombinant protein (Position: R33-R525). Human PGS1 shares 94.7% amino acid (aa) sequence identity with mouse PGS1.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PGS1 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.