

## Anti-RRP1B Antibody Picoband® Biotin Conjugated

Catalog Number: A08159-1-Biotin

### About RRP1B

RRP1B (ribosomal RNA processing 1 homolog B), also known as KIAA0179, is a human gene which is located on Chromosome 21. Enables transcription coactivator activity. Involved in several processes, including cellular response to virus; positive regulation by host of viral transcription; and positive regulation of transcription by RNA polymerase II. Located in chromosome; granular component; and nucleoplasm.

### Overview

Product Name	Anti-RRP1B Antibody Picoband® Biotin Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q14684

### Technical Details

Immunogen	E.coli-derived human RRP1B recombinant protein (Position: E20-K727).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-RRP1B Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.