

Anti-RRP1B Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A08159-1-Fluoro550

About RRP1B

RRP1B (ribosomal RNA processing 1 homolog B), also known as KIAA0179, is a human gene which is located on Chromosome 21. Enables transcription coactivator activity. Involved in several processes, including cellular response to virus; positive regulation by host of viral transcription; and positive regulation of transcription by RNA polymerase II. Located in chromosome; granular component; and nucleoplasm.

Overview

Product Name	Anti-RRP1B Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q14684

Technical Details

Immunogen	E.coli-derived human RRP1B recombinant protein (Position: E20-K727).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RRP1B Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.