

Anti-Frizzled 9/FZD9 Antibody Picoband® HRP Conjugated

Catalog Number: A08165-2-HRP

About FZD9

Members of the 'frizzled' gene family encode 7-transmembrane domain proteins that are receptors for Wnt signaling proteins. The FZD9 gene is located within the Williams syndrome common deletion region of chromosome 7, and heterozygous deletion of the FZD9 gene may contribute to the Williams syndrome phenotype. FZD9 is expressed predominantly in brain, testis, eye, skeletal muscle, and kidney.

Overview

| | |
|----------------------|---|
| Product Name | Anti-Frizzled 9/FZD9 Antibody Picoband® HRP Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. |
| Host | Rabbit |
| Uniprot ID | O00144 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived human Frizzled 9/FZD9 recombinant protein (Position: R27-L591). Human Frizzled 9/FZD9 shares 95.9% and 96.1% amino acid (aa) sequence identity with mouse and rat Frizzled 9/FZD9, respectively. |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | HRP |
| Suggested Dilutions | The intended application should be selected according to the customer's experimental requirements. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Frizzled 9/FZD9 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.