

Anti-Lymphocyte antigen 75 LY75 Antibody

Catalog Number: A08183-1


Overview

Product Name	Anti-Lymphocyte antigen 75 LY75 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Lymphocyte antigen 75 LY75 Antibody catalog # A08183-1. Tested in WB,IHC,IP applications. This antibody reacts with Human,Mouse.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O60449


Technical Details

Immunogen	Recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	IHC: 1:50-1:200 WB: 1:500

Anti-Lymphocyte antigen 75 LY75 Antibody (A08183-1) Images



Immunohistochemical analysis of paraffin-embedded human tonsil tissue using anti-LY75 antibody. Counter stained with hematoxylin.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue using anti-LY75 antibody. Counter stained with hematoxylin.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Lymphocyte antigen 75 LY75 Antibody

For Research Use Only. Not for use in diagnostic procedures.