

Anti-INSRR Antibody

Catalog Number: A08212-1


Overview

Product Name	Anti-INSRR Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-INSRR Antibody catalog # A08212-1. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P14616


Technical Details

Immunogen	Recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500 IHC: 1:50-1:200

Anti-INSRR Antibody (A08212-1) Images



Immunohistochemical analysis of paraffin-embedded human liver cancer tissue using anti-Insulin Receptor R antibody. Counter stained with hematoxylin.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue using anti-Insulin Receptor R antibody. Counter stained with hematoxylin.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-INSRR Antibody

For Research Use Only. Not for use in diagnostic procedures.