

Anti-SCAMP3 Antibody Picoband® HRP Conjugated

Catalog Number: A08215-1-HRP

About SCAMP3

Secretory carrier-associated membrane protein 3 is a protein that in humans is encoded by the SCAMP3 gene. This gene encodes an integral membrane protein that belongs to the secretory carrier membrane protein family. The encoded protein functions as a carrier to the cell surface in post-golgi recycling pathways. This protein is also involved in protein trafficking in endosomal pathways. Two transcript variants encoding different isoforms have been found for this gene.

Overview

Product Name	Anti-SCAMP3 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O14828

Technical Details

Immunogen	E.coli-derived human SCAMP3 recombinant protein (Position: M1-P347). Human SCAMP3 shares 90.3% amino acid (aa) sequence identity with mouse SCAMP3.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SCAMP3 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.