

## Anti-PSGR (F237) OR51E2 Antibody

Catalog Number: A08220


### Overview

Product Name	Anti-PSGR (F237) OR51E2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PSGR (F237) OR51E2 Antibody catalog # A08220. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9H255


### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 202-250 of Human PSGR.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

## Anti-PSGR (F237) OR51E2 Antibody (A08220) Images



Western blot (WB) analysis of PSGR (F237) polyclonal antibody at 1:500 dilution Lane1:Jurkat whole cell lysate Lane2:Raw264.7 whole cell lysate Lane3:PC12 whole cell lysate



Western blot (WB) analysis of PSGR (F237) pAb at 1:500 dilution Lane1:PC3 whole cell lysate(40ug) Lane2:L02 whole cell lysate(40ug) Lane3:AML-12 whole cell lysate(40ug) Lane4:H9C2 whole cell lysate(40ug)

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PSGR (F237) OR51E2 Antibody

For Research Use Only. Not for use in diagnostic procedures.