

Anti-LAS1L Antibody

Catalog Number: A08263-1

About LAS1L

LAS1L is a nucleolar protein required for the synthesis of the 60S ribosomal subunit and maturation of the 28S rRNA. LAS1L interacts with PELP1, TEX10, and WDR18, along with NOL9 and SENP3, to form a novel nucleolar complex that cofractionates with the 60S preribosomal subunit.

Overview

Product Name	Anti-LAS1L Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-LAS1L Antibody catalog # A08263-1. Tested in WB, IHC, ICC, IF, IP, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9Y4W2

Technical Details

Immunogen	E.coli-derived human LAS1L recombinant protein (Position: W3-T351).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400 ImmunoPrecipitation, 1:250-300 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LAS1L Antibody

For Research Use Only. Not for use in diagnostic procedures.