

Anti-Granzyme M/GZMM Antibody Biotin Conjugated

Catalog Number: A08286-2-Biotin

About GZMM

Human natural killer (NK) cells and activated lymphocytes express and store a distinct subset of neutral serine proteases together with proteoglycans and other immune effector molecules in large cytoplasmic granules. These serine proteases are collectively termed granzymes and include 4 distinct gene products: granzyme A, granzyme B, granzyme H, and the protein encoded by this gene, granzyme M. Two transcript variants encoding different isoforms have been found for this gene.

Overview

Product Name	Anti-Granzyme M/GZMM Antibody Biotin Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P51124

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human Granzyme M/GZMM. Human Granzyme M/GZMM shares 85.7% amino acid (aa) sequence identity with rat Granzyme M/GZMM.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Granzyme M/GZMM Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.