

## Anti-SLC26A8 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A08288-1-Fluoro647

### About SLC26A8

Predicted to enable chloride transmembrane transporter activity and organic anion transmembrane transporter activity. Predicted to be involved in chloride transmembrane transport; oxalate transport; and sulfate transmembrane transport. Predicted to be located in sperm annulus. Predicted to be active in plasma membrane. Is expressed in embryo. Human ortholog(s) of this gene implicated in spermatogenic failure 3. Orthologous to human SLC26A8 (solute carrier family 26 member 8).

### Overview

Product Name	Anti-SLC26A8 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, IF, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q8ROC3

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human SLC26A8.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-SLC26A8 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.