

Anti-Neuronatin/NNAT Antibody Biotin Conjugated

Catalog Number: A08311-1-Biotin

About NNAT

The protein encoded by this gene is a proteolipid that may be involved in the regulation of ion channels during brain development. The encoded protein may also play a role in forming and maintaining the structure of the nervous system. This gene is found within an intron of another gene, bladder cancer associated protein, but on the opposite strand. This gene is imprinted and is expressed only from the paternal allele.

Overview

Product Name	Anti-Neuronatin/NNAT Antibody Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q16517

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human Neuronatin/NNAT. Human Neuronatin/NNAT shares 95.8% and 91.7% amino acid (aa) sequence identity with mouse and rat Neuronatin/NNAT, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-Neuronatin/NNAT Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.