

Anti-GARNL1 (Phospho-T736) RALGAPA1 Antibody

Catalog Number: A08369T736

Overview

Product Name	Anti-GARNL1 (Phospho-T736) RALGAPA1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-GARNL1 (Phospho-T736) RALGAPA1 Antibody catalog # A08369T736. Tested in IHC, IF, WB applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q6GYQ0

Technical Details

Immunogen	Synthesized peptide derived from human GARNL1 (Phospho-T736)
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Suggested Dilutions	IHC 1:50-200 WB 1:500-2000 IF 1:50-200

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GARNL1 (Phospho-T736) RALGAPA1 Antibody

For Research Use Only. Not for use in diagnostic procedures.