

## Anti-CNOT4 Antibody

Catalog Number: A08405-1

### Overview

|                      |                                                                                                                                                                                                                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-CNOT4 Antibody                                                                                                                                                                                                                                                                                                                                                                                                        |
| Reactive Species     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                                                                                                                          |
| Description          | Boster Bio Anti-CNOT4 Antibody catalog # A08405-1. Tested in WB, IHC, ICC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.                                                                                                                                                                                                                                                                            |
| Application          | ELISA, IF, IHC, ICC, WB                                                                                                                                                                                                                                                                                                                                                                                                    |
| Clonality            | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Formulation          | 500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol<br>This antibody is supplied in a stabilized formulation.<br>Compatibility with conjugation reactions depends on the chemistry of the conjugation method used.<br>For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | 12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.                                                                                                                                                                                                                                                                                       |
| Host                 | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Uniprot ID           | O95628                                                                                                                                                                                                                                                                                                                                                                                                                     |

### Technical Details

|                     |                                                                                                                                     |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| Immunogen           | E.coli-derived human CNOT4 recombinant protein (Position: K275-T570).                                                               |
| Form                | Liquid                                                                                                                              |
| Concentration       | 500 ug/ml                                                                                                                           |
| Purification        | Immunogen affinity purified.                                                                                                        |
| Suggested Dilutions | Western blot, 1:500-2000<br>Immunohistochemistry, 1:50-400<br>Immunocytochemistry/Immunofluorescence, 1:50-400<br>ELISA, 1:100-1000 |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

### Anti-CNOT4 Antibody

For Research Use Only. Not for use in diagnostic procedures.