

Anti-CLPTM1 Antibody

Catalog Number: A08423-2

About CLPTM1

Cleft lip and palate transmembrane protein 1 (CLPTM1), a multipass transmembrane protein, was identified due to the association of a chromosomal translocation with a cleft lip and palate (PMID: 9828125). CLPTM1 regulates epileptic seizures by modulating GABAAR-mediated inhibitory synaptic transmission (PMID: 36519412). Highly expressed CLPTM1L may contribute to a poor prognosis and increase invasion, proliferation and migration of oral squamous cell carcinoma (PMID: 34586883).

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-CLPTM1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-CLPTM1 Antibody catalog # A08423-2. Tested in WB, IHC, IP, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IP, IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O96005

Technical Details

Immunogen	E.coli-derived human CLPTM1 recombinant protein (Position: A3-A348).
Form	Liquid

Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 ImmunoPrecipitation, 1:250-300 ELISA, 1:100-1000

[Anti-CLPTM1 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.