

## Anti-Myosin light chain 4 MYL4 Antibody

Catalog Number: A08496


### Overview

Product Name	Anti-Myosin light chain 4 MYL4 Antibody
Reactive Species	Rat
Description	Boster Bio Anti-Myosin light chain 4 MYL4 Antibody catalog # A08496. Tested in WB,IHC applications. This antibody reacts with Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P12829

### Technical Details

Immunogen	Recombinant Protein of human MYL4
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:2000-1:5000 IHC: 1:50-1:200

## Anti-Myosin light chain 4 MYL4 Antibody (A08496) Images



Western blot analysis of extracts of rat heart, using MYL4 antibody at 1:3000 dilution.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Myosin light chain 4 MYL4 Antibody

For Research Use Only. Not for use in diagnostic procedures.