

Anti-Apolipoprotein C-IV APOC4 Antibody

Catalog Number: A08520


Overview

Product Name	Anti-Apolipoprotein C-IV APOC4 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Apolipoprotein C-IV APOC4 Antibody catalog # A08520. Tested in WB,IHC applications. This antibody reacts with Human,Mouse.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P55056


Technical Details

Immunogen	A synthetic peptide corresponding to residues in Human APOC4.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 105% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

Anti-Apolipoprotein C-IV APOC4 Antibody (A08520) Images



Western blot (WB) analysis of APOC4 polyclonal antibody at 1:500 dilution Lane1:HEK293T whole cell lysate Lane2:sp2/0 whole cell lysate Lane3:H9C12 whole cell lysate



Immunohistochemistry (IHC) analyzes of APOC4 pAb in paraffin-embedded human liver carcinoma tissue at 1:50, showing secreted staining. Negative control (the right) Using PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG-biotin followed by avidin-peroxidase.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Apolipoprotein C-IV APOC4 Antibody

For Research Use Only. Not for use in diagnostic procedures.