

Anti-Keratin, type II cuticular Hb6 KRT86 Antibody

Catalog Number: A08530

Overview

Product Name	Anti-Keratin, type II cuticular Hb6 KRT86 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Keratin, type II cuticular Hb6 KRT86 Antibody catalog # A08530. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O43790

Technical Details

Immunogen	Synthesized peptide derived from part region of human protein AA range: 416-466
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	KRT86 antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-2000 ELISA 1:5000-20000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Keratin, type II cuticular Hb6 KRT86 Antibody

For Research Use Only. Not for use in diagnostic procedures.